

1 A Graphical User Interface for Individualized Locomotor
2 Training of Infants With or at High Risk of Cerebral Palsy
3 Using a Robotic Assistive Device

4 Monique O. Shotande*¹, Julie Skorup², Stephanie O’Leary², Morgan Alcott²,
5 Beth A. Smith^{3,4,5}, Laura A. Prosser^{6,7}, Mustafa Ghazi⁸,
6 Thubi H. A. Kolobe^{7,9}, and Andrew H. Fagg^{1,10}

7 ¹School of Computer Science, University of Oklahoma, Norman, USA

8 ²Department of Physical Therapy, The Children’s Hospital of Philadelphia, Philadelphia,
9 USA

10 ³Developmental Neuroscience and Neurogenetics Program, The Saban Research Institute,
11 Los Angeles, USA

12 ⁴Division of Developmental-Behavioral Pediatrics, Children’s Hospital Los Angeles, Los
13 Angeles, USA

14 ⁵Department of Pediatrics, Keck School of Medicine, University of Southern California,
15 Los Angeles, USA

16 ⁶Department of Pediatrics, University of Pennsylvania, Philadelphia, USA

17 ⁷Division of Rehabilitation Medicine, The Children’s Hospital of Philadelphia,
18 Philadelphia, USA

19 ⁸Department of Pediatrics, Children’s Hospital Los Angeles, Los Angeles, USA

20 ⁹Department of Rehabilitation Sciences, University of Oklahoma Health Sciences,
21 Oklahoma City, USA

22 ¹⁰Institute for Biomedical Engineering, Science and Technology, University of Oklahoma,
23 Norman, USA

24 August 8, 2024

*Corresponding author: mo.shotande@gmail.com

Abstract

Background

Precision rehabilitation is essential to optimally meet the specific needs of children with developmental motor disorders such as cerebral palsy (CP). While assistive devices are useful tools for motor training by supporting the child, augmenting actions, and encouraging practice, these devices often have limited sets of predetermined, [fixed modes of interaction and parameters](#) that can limit the effect of their responses to performance.

Objective

We present [the design of a graphical user interface \(GUI\) that facilitates prone locomotor learning and monitoring of individualized learning efforts by infants using a robotic assistive crawling device.](#)

Methods

The Self-Initiated Prone Progression Crawler (SIPPC) robot engages sensor-driven interaction and incorporates motor learning and precision rehabilitation principles to encourage and assist early prone locomotion. [Training on the SIPPC is 12–16 weeks, with hands and feet contacting the floor. Infant and robot data are displayed in the GUI, where therapists monitor infant performance and robot behavior. Through the GUI, therapists can make real-time adjustments in the robot response to infant movements, thereby increasing opportunities for successful movement at any learning stage. Selectable robot control modes either respond to infant-generated ground reaction forces \(force control\) or limb gestures \(G1 and G2\). The GUI enables control mode response threshold and magnitude adjustments, and provides scripted session guidance and data management.](#)

Results

[Preliminary results show therapists can provide real-time configuration changes to the SIPPC response type, sensitivity, and magnitude based on infant skill, mood, and effort. Infant and robot behavior are observed concurrently with these configurations. Therapists enable force control and G1 for most sessions. G2 detects more coordinated limb gestures than G1 and used during later sessions after infants have more practice and SIPPC experience. During early therapy sessions, therapists change parameters multiple times throughout the trials. In subsequent weeks, therapists typically adjust parameters at the trial beginning.](#)

Conclusion

Individualized intervention is important for optimizing rehabilitation outcomes. Our interface [design allows therapists to customize the assistive robotic responses to each infant's needs and learning style to promote motor learning specific to each infant's motor function. This innovative approach allows precise monitoring of movement learning across time and infants, offering additional value to measuring infant locomotor development.](#)

Keywords

motor learning, motor training, motor control, assistive therapy robotic systems/devices, crawling skill acquisition, prone progression, cerebral palsy (CP)

61 Introduction

62 Cerebral Palsy (CP) is a group of developmental motor disorders in infants that results in lifelong impairment
63 and can impact the progression of locomotion, cognition, perception, communication, and behavior [1–3]. Impaired
64 motor control and coordination in infants with CP can vary across individuals and produce challenges
65 to motor learning and skill acquisition, leading to delayed or missed developmental milestones, such as crawling
66 ing [4, 5]. Early diagnosis and intervention are therefore crucial for infants at high risk of CP to leverage neural
67 plasticity and neural pathway pruning early in life [2–7].

68 Crawling is a profound milestone in development for most typically developing infants; acquisition of this
69 milestone influences muscular strength, motor coordination, and cognitive development [8–10]. Typically,
70 infants learn to crawl by moving through a sequence of developmental stages, progressively increasing the
71 magnitude and coordination of their movements [11, 12]. When placed prone, they initially learn to lift their
72 heads for increasing duration to observe their surroundings and to identify objects of interest. As they obtain
73 more strength and control over their trunk and limbs, they start to shift their weight and reach for nearby
74 objects. With increasing muscle strength and coordination, they start moving around on their stomachs, before
75 transitioning to [quadrupedal motion](#) on their hands and knees. With further refinement of limb and trunk control
76 and coordination, they can move in their environment toward toys and other objects of interest. [However,](#)
77 [challenges such as reduced muscle strength or control can result from disorders like CP and make achieving](#)
78 [any of these stages unassisted much more challenging \[1–3, 13\].](#)

79 [Activity-based interventions can potentially improve motor function by targeting multiple motor perfor-](#)
80 [mance components, such as endurance, coordination, and strength, with varied levels of intensity and frequency](#)
81 [to provide individualized rehabilitation \[3, 4, 7, 14\]. For example, passive range motion exercises can help](#)
82 [prevent or reduce joint contractures and low resistance high repetition exercises can help with muscle strength-](#)
83 [ening and endurance \[3\]. Morgan et al. \[7\] developed the Goals-Activity-Motor Enrichment \(GAME\) therapy](#)
84 [approach to promote early, intense, environmentally enriched, task-specific, home-based training for motor](#)
85 [and cognitive skill improvement in infants with CP. In the GAME paradigm, therapists and caregivers work](#)
86 [together to adjust the infant’s environment to promote successful self-generated motor activity. While tem-](#)
87 [porary performance improvements are motivating in therapy, it is crucial to ensure advancements in motor](#)
88 [capabilities are permanent and applicable to active daily life. To do so, is useful to increase motor challenge as](#)
89 [performance improves to help solidify learned problem-solving capabilities and skill generalization \[7, 15, 16\].](#)

90 [While therapeutic interventions can aid individuals in learning motor skills to navigate their environments](#)
91 [or to learn tasks necessary for daily activities \[17, 18\], there is often a high variability in individual responses](#)
92 [to intervention \[19\] and some variation in therapist implementation of intervention procedures \[13, 14\].](#)
93 [Therefore, interventions should consider the unique dispositions of each individual and be adjustable to best](#)
94 [address their personal needs and limitations. And, these adjustments should be quantitative to allow for](#)
95 [consistency and objective comparison between alternative therapist implementations.](#)

96 [Figure 1 summarizes a model of the dynamic, real-time interaction between the environment, device,](#)
97 [learner, and therapist. Through a feedback loop, the device augments the learner’s capacity to explore](#)
98 [the space of relevant constituent motor actions with differing complexities. The therapist can observe the](#)
99 [activities of the learner and respond by altering the learner’s environment. In addition, the device provides](#)
100 [detailed information about the performance of the learner to the therapist, who, in turn, tunes the response](#)
101 [of the device to meet the learner’s specific needs.](#)

102 [In this context, an assistive therapy device that applies to a range of motor deficits of varying severity](#)
103 [levels should adhere to the following *design principles*:](#)

- 104 1. support learning of a range of skills of varying complexities,




Figure 1: Flow of extrinsic feedback. In natural environments, there is a feedback loop between the environment and the learner. Without an assistive therapy device, the therapist acts as a mediator between the environment and the learner by controlling the configuration of the environment and how the learner interacts with it. With a therapy device, both the learner and the therapist obtain feedback from the device, and the therapist can modify the configuration of that feedback. The assistive therapy device acts as an additional mediator with the environment. Additionally, the device can respond in real-time to the performance of the learner without waiting on input from the therapist. In this study, the SIPPC is the device and the infant is the learner.

- 105 2. provide automatic extrinsic feedback to the learner through repeatable assistance that is contingent on
 106 the learner's performance in real-time, and
- 107 3. provide extrinsic feedback to the therapist to allow objective performance assessment and to adjust the
 108 learner-device interaction.

109 The first design principle addresses the variation in motor capacity of individuals with motor disorders. As-
 110 sistive devices can augment and provide some enhancements to individual movements. This can help learners
 111 practice and build new motor skills. Additionally, the device helps the learner identify and engage effective
 112 skills related to desired motor objectives. Because CP can impact motor capability in many ways, it is crucial
 113 to have devices that are capable of supporting the practice and achievement of a large set of skills appropriate
 114 to the individual's level of motor function.

115 The second design principle requires that the device provide extrinsic feedback to the learner in real-time
 116 in response to the learner's performance. Moinuddin et al. [20] defines feedback, in the context of motor
 117 learning, as information "fed back" to the learner as the result of movement. The feedback can be sent before,
 118 during (real-time) or after the action to inform motor adjustments in preparation for subsequent actions. Real-
 119 time feedback is important to allow the learner to quickly change strategy and to determine the relationship
 120 between certain actions and progress towards a goal [4, 16]. Solidification of learning through quality repetitive
 121 practice is a core component of effective therapy and motor learning [17, 18, 20]. Salmoni et al. [15]
 122 and Kitago and Krakauer [16] discuss the temporary impacts of extrinsic feedback (knowledge of results;
 123 KR), and that reducing or removing KR in later stages attenuates the temporary performance effects of
 124 the feedback and accentuates lasting, transferable improvements that are actually due to learning. This
 125 form of real-time extrinsic feedback can provide infants at risk of CP with high-frequency training and help
 126 implicitly communicate learning objectives to infants since they are not able to follow instructions and verbally

127 communicate.

128 The third design principle requires that the device provides qualitative and quantitative feedback to the
129 therapist. It is paramount that therapists have objective knowledge about the status, progress, and trajectory of
130 the learner. This knowledge empowers the therapists to make informed decisions for the learner's experience.
131 Upon receipt of the feedback, the therapists can alter the feedback response to future actions during the same
132 session. Modifying the feedback in real-time enhances the therapist's ability to affect behavior and guide
133 learning. Changes in feedback can serve to motivate activity or encourage the exploration of various sets of
134 underutilized actions (e.g., encouraging infants who struggle to move one particular arm to begin engaging it
135 to reach for toys). As the learner improves, the therapist should be able to modify the feedback to increase
136 the challenge and prevent primary dependency on the feedback. Moreover, if the learner is not successful,
137 the therapist can decrease the challenge to provide successful feedback to the learner.

138 Assistive devices (e.g., walkers, gait trainers, exoskeletons, robotic devices) can be used in therapy to
139 augment and enhance motor performance by enabling exploration and discovery in various environments [3,
140 13, 14]. The combined extrinsic feedback from the environment and an assistive device can help learners
141 adapt, evaluate the utility of various motor strategies [21], and improve overall motor function [3]. Assistive
142 therapy devices can also be useful tools for augmenting the interaction between the therapist and the learner.
143 Through a combination of sensors, functional user interfaces, and robotics, these devices have the potential
144 to offer therapists more quantitative information, control, and flexibility when making therapy decisions, while
145 simultaneously providing feedback and assistance to train the learner. However, such devices are often limited
146 in applicability and availability for children, especially infants, due to differences in needs between children
147 and adults. Moreover, assistive devices for infants that support early learning of locomotor skills, such as
148 crawling, are not widely available [22, 23], do not provide performance feedback [24, 25], or do not offer a
149 flexible set of parameterized controls for the therapist [24–26].

150 The Self-Initiated Prone Progression Crawler (SIPPC) robot is a unique assistive crawling device that is
151 designed to promote movement learning in infants with motor delays by offering an array of crawling-based
152 movement strategy options [21–23, 27]. The SIPPC supports the weight of the infant in a prone posture in
153 which their hands, feet and knees can contact the floor. The robot can carry the infant along a combination
154 of forward and turning paths – contingent on the ground reaction forces (GRFs) or the limb movements (LMs)
155 produced by the infant. For example, if the infant pushed against the floor to the left with their two hands,
156 the SIPPC would respond by turning the infant to the right (GRF). Or, if the infant made a pulling motion
157 with their hands from in front of their body toward themselves, even without contacting the floor, the SIPPC
158 would move the infant forward (LM).

159 Kolobe and Fagg [21] investigated the effects of these response contingencies with three groups of infants:
160 typically developing infants with GRF-based responses (TD), infants at risk for CP with GRF-based responses
161 (GRF), and infants at risk for CP with both GRF and LM responses (GRF+LM). Over bi-weekly exposure
162 to the SIPPC for 12–20 weeks, the TD group demonstrated a greater increase in the magnitude of arm
163 movements and in trial-and-error activities over both at risk groups. In addition, among the at risk groups,
164 the GRF+LM group demonstrated an advantage in these metrics over the GRF group. This indicates that
165 the SIPPC with appropriately selected response contingencies can contribute to motor exploration and to
166 learning new motor skills. However, the fixed contingencies by group can present limitations on an individual
167 basis. Adding the LM contingency on top of the GRF contingency can make the SIPPC easier for the infant
168 first to learn to engage the robot in producing purposeful movements, but does not necessarily challenge the
169 infant to refine their skills further. On the other hand, the GRF-only contingency is more challenging to the
170 infant, but to a degree that can impede the first learning steps.

171 For the SIPPC (and other therapeutic devices) to meet the unique needs of more individuals, it is important
172 to provide the therapists with the ability to dynamically select contingencies, in addition to tuning their details,
173 in an individual and context-dependent manner. A key challenge is choosing this set of contingencies and

174 their associated parameters, and then presenting these options to the therapist in a manner that facilitates
175 (and not impedes) their interaction with the patient. Many existing therapy devices possess user interfaces
176 that operate as a control center for the therapy session to adjust the exercises based on therapy goals, the
177 patient needs [28–30], and they take steps to promote understandability, reliability, and accessibility [30, 31].
178 Octavia et al. [28] designed Matti, an interactive gaming mat that uses a tangible user interface (TUI) to
179 support motor rehabilitation for children with developmental coordination disorder (DCD). Here the TUI is
180 designed to engage and motivate the patient and a screen displays visual feedback from therapy games. A
181 separate graphical user interface (GUI) is provided to the therapists to select games and adjust configurations
182 based on the therapy goals. One limitation of this design is that once the therapist begins the therapy game,
183 the configuration is fixed until the game is complete or terminated. This can limit practice opportunities
184 that optimally address the child’s unique motor function and skill level. In work by Tucker et al. [29], they
185 developed a GUI for wireless operation and data logging from a robotic knee exoskeleton for children with CP.
186 The interface allows users to select and configure two control modes and automatically enforces the order of
187 procedural therapy operations to reduce the chance of error. Despite the capabilities of the exoskeleton and
188 GUI, the accessibility of the device in Tucker et al. [29] is limited since access to the GUI requires an intricate
189 installation procedure.

190 In this paper, we describe the design of the SIPPC’s graphical user interface (GUI) and its application of
191 these design principles. Through the GUI, we can offer a fine-level of control of the therapy device such that
192 the device can be used to support a larger set of patient needs and facilitate precision rehabilitation. We also
193 present preliminary data from an ongoing study of infant locomotor learning.

194 **Materials and System**

195 **Self-Initiated Prone Progression Crawler**

196 The SIPPC 4 (Figure 2) is the latest generation in a series of integrated assistive crawling robotic systems that
197 promotes prone locomotor learning in infants with motor delays [21]. The infant lies prone on a padded platform
198 such that their hips, trunk and potentially their head are supported, while their hands, feet and knees are able
199 to contact the floor. The SIPPC uses information from a set of sensors to assist the infant in motor training
200 for precision rehabilitation. *Embedded inside the platform is a 6-degrees-of-freedom force/torque (F/T) sensor*
201 *with which ground reaction forces and weight shift can be measured.* An infant motion capture suit tracks limb
202 and trunk movements in real-time using twelve Inertial Measurement Units (IMUs) [32]. The 3D orientation
203 data from the IMUs are combined with a skeletal model of the infant to estimate the 3D locations of key points
204 on the infant’s body, including the wrists, ankles, and feet [33]. Three cameras capture lateral views of the
205 infant, the infant’s head, hands and feet, as well as area in front of the infant. The images from the cameras are
206 composited into a single view by a dedicated video processor; the composite image is recorded to disk and sent
207 to the connected client devices (e.g., tablet, smartphone).

208 The SIPPC has two drive wheels mounted on either side of the infant’s hips. Each wheel is independently
209 controlled by a local microcontroller that communicates with the control processor, which is responsible for
210 coordinating their behavior and monitoring their state for safety. A rear ball castor wheel sits behind the infant
211 as a third point of contact and a front skid plate prevents the SIPPC from tipping. Flexible platform covers at
212 the front of the platform and around the wheel pods prevent the infant’s hands from being pinched when the
213 robot moves. *Additional details about the hardware design of the SIPPC are described by Ghazi et al. [34].*
214

215 **Software System Design** The SIPPC software system is built on top of the Robot Operating System (ROS,




Figure 2: A therapy session where the infant is using the Self-Initiated Prone Progression Crawler (SIPPC 4). The parent is in front of the infant encouraging her to move her limbs and progress forward. The infant is outfitted with an eleven-segment motion capture suit. Three cameras record the limb movements and the space in front of the infant. Under the red platform, a force/torque sensor measures the ground reaction forces (GRFs) that she generates. The yellow box behind the infant is the electronics box housing most of the robot's processors, the battery, and the power regulation system.

216 v1, Kinetic Kame) [35], which supports the flexible, real-time communication and coordination between a
 217 large number of sensory, logging and control processes (referred to as ROS nodes). These ROS nodes are
 218 distributed across multiple processors and implemented with a mix of C++11 and Python. The base update
 219 rate for the real-time processes (including sensing, control and logging processes) is 50Hz. Event-based
 220 processes wait for a trigger event to perform an update, such as parameter change events from the GUI.

221 The software architecture details are outlined in Figure 4 and described by Ghazi et al. [34]. In short,
 222 the sensor nodes gather and publish data from the force/torque sensor, the motion capture suit, the SIPPC
 223 wheel state, the cameras, and the battery. A low-level control node is responsible for producing smooth
 224 and safe motions of the SIPPC. A video compositor node is responsible for the real-time integration of the
 225 images from three separate cameras and graphical representations of the SIPPC state into a single video
 226 stream, which is then transmitted to the GUI at 10FPS and recorded to disk for later coding. A session
 227 control node is responsible for enforcing the sequence of allowable steps that the therapist can take during a
 228 recording session. A logging node is responsible for writing all key SIPPC state information to a log file for
 229 later analysis. And, a processor node is responsible for monitoring the CPU load, memory, and data storage
 230 on each processor.

231 The GUI communicates through ROS to present the robot system state to the therapist and to enable
 232 the therapist to tune the robot responses to each infant's personal needs. The GUI receives regular robot and
 233 infant state updates through ROS at 10Hz. The therapist is also provided with diagnostic details for real-time
 234 troubleshooting should issues arise. ROSBridge [36] is a ROS package enabling communication between the
 235 ROS environment on the robot and the client devices (e.g., phone, tablet). The package uses WebSockets to
 236 transmit data in the JavaScript Object Notation (JSON) [37] format and maintain a shared robot state across
 237 all connected devices. The robot uses a RaspberryPi (RPi) [38] as an access point (AP) for client devices to




Figure 3: Web Architecture. Robot data flows between the ROS environment on the robot (using WebSockets) and client devices (e.g., phone, laptop, tablet). The Flask web server on the robot hosts and transmits the pages to the connected clients. State change requests flow between the Flask server on the robot and the client device and back to ROS to valid state change requests. The ROSBridge server carries the protocol for converting ROS messages into JavaScript Object Notation (JSON) for transport through the WebSocket.

238 connect to and a web server for users to access the various pages of the user interface. A high-level diagram
 239 of the web architecture is shown in Figure 3. The web application server uses Python Flask [39] to define the
 240 endpoints corresponding to each of the interface's pages. Each page provides information about the infant's
 241 progress, controls for monitoring the robot's state, or capabilities for managing the data. Guicorn [40] is a
 242 gateway server that manages multiple processes for handling client requests in parallel when multiple users
 243 are simultaneously connected to the interface. Nginx [41] acts as a reverse proxy server for port forwarding.
 244 For dynamic device layout adjustment, Bootstrap generates responsive styles that allow simple specifications
 245 for automatic layout rearrangement based on the size or resolution of the client's device.
 246

247 **Networking** An outline of the network and device connections is shown in Figure 5. The robot has
 248 five specialized processors (three RaspberryPi's (RPi) and two Mbed microprocessors). These processors
 249 communicate on the same subnet via Ethernet. The main RPi (sippcbridge) acts as an access point and web
 250 server. Client devices (e.g., tablet) connect to sippcbridge to access the SIPPCC's local password-protected
 251 WiFi network and the GUI. The main RPi also hosts the ROS master process, which manages the ROS nodes
 252 and the inter-node communication. The robot can be connected to the internet through the sippcbridge to
 253 synchronize data between the secure cloud storage, other SIPPCC robots, and secure remote lab devices for
 254 software updates and maintenance. The video RPi (sippcvideo) is dedicated to capturing video from the
 255 cameras and compositing them into one video stream along with data plots of the kinematic and robot data.
 256 The solid state drive (SSD) for storing all the session data on the robot is connected to sippcvideo. The SSD
 257 is mounted by the other RPi's over the robot's local network. The control Mbed (sippccontrol) monitors and
 258 manipulates the wheels and force/torque (F/T) sensor. The suit Mbed is connected to the twelve IMUs and
 259 is responsible for collating and transmitting the 3D orientation data to other nodes in the network.

260 **Robot Assistance Control Modes**

261 Control modes are types of robot assistance that trigger responses by the robot based on actions that the infant
262 performs. There are two general types of control modes the therapists can engage for the robot: proportional
263 and discrete. For the force control mode, the robot acceleration is linearly related to the magnitude of the
264 infant-generated GRFs. The gain of this relationship is a parameter that the therapist can adjust under the
265 Force/Torque Control panel on the Parameter Page (Figure 7; Forward Force Boost and Turning Torque
266 Boost parameters).

267 The discrete control modes produce temporally extended, minimum jerk trajectory profiles in the forward,
268 left or right directions. Only one discrete action can be triggered at a time; its acceleration profile is layered
269 on top of that of the force controller. The parameters of the profile (temporal duration and distance) can
270 be altered by the therapist. Several different discrete control modes are available: power steering (PS), suit
271 gesture recognition, and manual control. Assistance from PS is triggered when the GRFs exceed a threshold
272 that the therapist defines. Two suit gesture control modes are available: Gesture 1 (G1) and Gesture 2 (G2).
273 G1 is a heuristically defined set of rules that translate limb motion into triggering a discrete robot response [32].
274 G2 is a separate set of machine-learned (ML) gesture rules encoded with a forest of relational probability
275 trees [42, 43]. This model was trained using kinematic and PS assistance event data from typically developing
276 infants learning to use the SIPPC with PS only [21]. G2 identifies additional, more precise types of limb
277 movements beyond what can visually be observed and defined by the therapists. G1 responds more to coarse-
278 level movements in the arms and legs, whereas the sensitivity of G2 can be dynamically adjusted to encourage
279 more coordinated crawling-like movements. Specifically, the G2 model produces as output an estimated
280 probability of a forward, left or right movement. When one of these predictions exceeds a therapist-defined
281 threshold, the G2 control mode generates the corresponding discrete movement. As a probability threshold
282 increases, the infant must produce movements that better match the movements of the typically developing
283 infants from the training set. These parameters are under Session Configuration: Gesture 2 Sensitivity on
284 the Session View page (Figure 6). Finally, the Manual control mode allows the therapist to remotely drive the
285 robot.

286 **User Interface**

287 The SIPPC GUI is used by the therapist to guide motor training by allowing the therapists to control the difficulty
288 of the sessions, dynamically. The GUI allows the therapist to adjust the control parameters of the robot in real-
289 time to modify its responsiveness based on the infant's engagement and movement. The therapist can modify
290 multiple parameters as the infant becomes more skilled at driving the robot to ensure that the infant is continually
291 challenged to practice and learn more difficult motor skills. Therapists can adjust parameters such as:

- 292 • GRF response threshold to adjust performance difficulty,
- 293 • limb and trunk gesture recognition sensitivity to adjust performance difficulty,
- 294 • robot control modes to target limb use for generating movement, and
- 295 • robot movement magnitude and duration to motivate limb and trunk movement.

296 Early on, when the infants are less skilled, the robot is configured to be very responsive. Sessions can be made
297 less challenging for infants with more severe motor impairments or more encouraging for infants with less
298 motivation to move. This allows infants to practice motor activities at their current skill level. For example,
299 if the infant is visually engaged with a toy but not moving towards it, the therapist can make it easier for the
300 infant to achieve success. The therapist can increase the gesture sensitivity to detect small limb movements

301 of the infant to drive the robot toward the toy. Conversely, if the infant is easily generating robot movements
302 by actively reaching, pushing, kicking, and progressing toward toys, the therapist can increase the challenge
303 of driving the robot. This can be achieved with multiple combinations of changes to the robot parameters in
304 the GUI. For example, the therapist can decrease the distance that the robot travels and increase the amount of
305 force that the infant must generate to move toward the desired item.

306 The SIPPC GUI promotes the proposed design principles by being accessible by most personal devices (e.g.,
307 tablet, laptop, smartphone, etc.), maintaining a consistent shared state across all connected devices, being self-
308 contained and informative, being modular in design, and limiting session procedural operations to a specified
309 sequence. To meet these requirements, the GUI was built as a local web application. Web applications are
310 accessible on any device with WiFi connection abilities and are hardware agnostic, eliminating the need to
311 explicitly develop the interface for multiple types of devices. Additionally, web applications do not require
312 installation on a client device and are accessed simply by connecting to the robot's local network and the
313 SIPPC URL from a browser. For security, the robot and its local WiFi network are password protected. One
314 must connect to the robot's network to access the GUI.

315 Most operations in the GUI are performed by pressing a single toggle button and changing parameters
316 is done by selecting checkboxes or moving a slider. Additionally, the interactive components of the interface
317 (e.g., buttons) change appearance if interacted with correctly, to clearly indicate to the users the button's
318 current state. This simplifies the interface and the workflow for the therapists. For example, the toggle
319 buttons change color and text to indicate state information. In particular, the button to start and stop Data
320 Recording is red and displays "STOPPED" when data recording is not running. When data recording is
321 running, the color changes to yellow and the text inside the button is set to "STARTED." The sliders display
322 the value for their corresponding parameter in a text box that the user can use to type the value instead if
323 typing is preferred to using the slider.

324 Multiple client devices can be connected to the GUI of the same robot simultaneously to observe the system
325 state. To ensure a consistent state across the connected devices, each device receives the state data from a
326 centralized parameter ROS node; any parameter changes from the therapist are sent as change requests to this
327 node, which validates the change before altering the system state. The ability to have multiple client devices
328 connected to the robot's interface also assists therapists with training others to use the system, involving the
329 infant's caregivers and monitoring safety information. The user interface is separated into eight pages (Session
330 View, Parameter Page, Post Processing, Subject Page, Clinician Page, User Guide, System Details, and Admin
331 Dashboard; Figure 6 and 7). Each page focuses on specific therapist tasks. Troubleshooting details and a user
332 guide are included in the interface for convenience to clarify how various interface tools operate and to explain
333 the causes and solutions to issues that might arise while users are interacting with the SIPPC.

334 **Session View** The main interface page is the Session View (Figure 6). This page displays the robot and sensor
335 status information (e.g., whether the force sensors are connected), and allows the therapist to engage/disengage
336 the robot, drive the robot (i.e., Manual Assist panel in Figure 6), perform suit calibration for the trunk and limbs
337 (Suit Calibration panel in Figure 6), modify the session configuration (i.e., robot-infant interaction parameters),
338 and control the trial state (e.g., start or pause the trial timer). The therapist can also take notes that are logged
339 along with other system data for use during subsequent sessions and analyses.

340 The panels on the page also facilitate modularity in the interface design. The panels are automatically
341 rearranged for different screen resolutions and are laid out in an order that reflects the typical temporal order
342 of the activities taken by a therapist during a session/trial. The interface will not allow certain actions to be
343 performed before others. If the user requests an invalid sequence the interface notifies the user of the issue
344 and the correct action to perform. For example, before Data Recording can start, the connection to the wheels,
345 force/torque sensors, and other processors must be established (System Status panel in Figure 6, each indicated
346 by a green check mark), and the subject and clinician must be selected. Precautionary checks are built into the

347 SIPPC. The robot cannot be engaged without the wheels and force sensors being connected. Before starting the
348 first trial, the suit must be connected and the suit and the Force/Torque Sensor must be calibrated. Therapists
349 can start the trial and engage the robot and the infant's control when all precautions are met. This design reduces
350 errors in the procedures and promotes consistency across sessions, subjects, and therapists. The checks in the
351 system also reduce the chances of erroneous interactions with the robot that could result in damaged or missing
352 data.

353 **Parameters Configuration** The Parameters Page (Figure 7) exposes robot-infant interaction parameters, al-
354 lowing in-session adjustments. Modifications to these parameters enable the therapists to guide motor training
355 on an individual basis. These configurations can be saved and shared across multiple robots on an individ-
356 ual basis. The Discrete Assistance panel modifies the parameters impacting the robot's response to assistance
357 events triggered by the power steering or gesture recognition control modes. The Power Steering panel allows
358 therapists to adjust the threshold GRFs necessary to provide the infant with extra force-based assistance to
359 move forward, left, or right. The lower the threshold, the easier it is to trigger a power-steering-assisted boost.
360 Force/Torque Sensor parameters manage the filtering of sensor noise. The Force/Torque Control parameters
361 handle the robot's level of continuous response to the infant's GRFs (the Global Gain and Boost parameters).
362 Increasing the Boost parameters increases the robot's acceleration proportional to the GRFs and can provide
363 infants with extra motivation in response to smaller forces when the infant might be struggling to produce
364 forces. Lastly, the damping parameters control how quickly the robot slows down after a movement has been
365 completed.

366 **Data Handling** The Post-Processing page (Figure S1), guides users through the steps for data post-processing,
367 uploading, and deletion. The structure and functionality of the page restrict the order in which the steps can be
368 performed to prevent loss of data and minimize the time and effort spent on data management. For example,
369 video data cannot be uploaded until after it is compressed. Files cannot be deleted until after they have been
370 uploaded. Therapists can also download subject parameter configuration files and notes from other robots and
371 sessions to use for current sessions, allowing the different robots to be interchangeably used for therapy ses-
372 sions. All data are uploaded to the HIPAA-compliant version of Box. Only authorized users can access this
373 data repository. All data exchanges with Box are encrypted in transit, and Box encrypts the data at rest.
374 Post-processing is done on a computing system at the University of Oklahoma (OU) where the data are stored
375 on a ZFS-based system with encrypted disks and three parity disks to prevent data loss. Data exchanges
376 between this file system and the processing computers are also encrypted in transit. Data (including video)
377 are stored on the SIPPC until the videos can be transmitted, processed, and human-verified (typically takes
378 one week). See the Supplemental Details for details on all the pages in the GUI.

379 **Methods**

380

381 **Design Requirements**

382 The current SIPPC robot and interface evolved from prior studies involving the SIPPC 2 [21, 22, 44] and
383 SIPPC 3 [23, 27]. SIPPC 2 was the first functional prototype to be used in a study that allowed for
384 the comparison of typically developing and at risk infants. The robot passively accelerated in response to
385 GRFs, but the control system could also generate short duration acceleration "bumps" in response to small
386 passively-generated movements (similar to PS) and to gesture recognition events (G1, specifically). The

Goal	Requirement	GUI Page
Facilitate personalization of therapy and monitoring	Select infant and clinician profiles	Session View
	Start/stop data recording sessions	
	Start/pause/resume/stop trials within each data collection session	
	Calibrate the suit and force/torque sensors	
	View infant performance and robot behaviors	
	Control robot response and behavior (i.e., ability to change control parameters)	
Facilitate monitoring	View robot system state and details (e.g., battery, sensor states, processor load and memory, data storage space, etc)	Session View and System Details
	Provide real-time troubleshooting information and robot and GUI documentation	All and User Guide
	Enforce compliance with the study procedures and protocol	All Pages
	Accessible from network-enabled devices (e.g., phone, tablet, laptop, etc)	
Data Management	Share infant profile data between robots such that they are interchangeable	Post-Processing
	Manage data (upload, download, delete)	
	Shut down robot	System Details

Table 1: GUI Design Requirements

387 wire-connected interface computer provided the therapist with minimal abilities to start/stop trial recording
388 to visualize the state of a small battery of sensors. SIPPC 3 expanded the sensory and motor capabilities
389 of the system, allowing for higher fidelity data collection and the ability to prototype additional movement
390 contingencies, including PS, G1 and G2. The GUI executed on a wireless laptop computer, and allowed
391 the therapist to select control modes and adjust a small number of control parameters. However, following
392 the prototype phase, these parameters were fixed for the purposes of the subsequent study that compared
393 typically developing infants with at risk infants.

394 SIPPC 4, the subject of this paper, is part of a larger study that is focused on providing personalized motor
395 training to the infants using fine-level dosage control [45]. The new requirements for this therapy-focused
396 study include: a flexible interface that is accessible from any WiFi-enabled device, and GUI-based real-time,
397 fine-level control over the device responses, as well as the data management processes. A key challenge in
398 allowing the therapists to have fine-level control of the SIPPC behavior was in balancing this complexity with
399 their need to focus on the infant during the therapy sessions. Table 1 summarizes the specific requirements
400 to achieve these goals.

401

402 **Validation and Testing**

403 Prior versions of the SIPPC allowed us to prototype a range of ideas for supporting infant crawling learning and
404 infant-robot interaction. The SIPPC 4 represents a consolidation of these efforts into a device that is practical
405 for use in homes and flexible enough to allow for a wide range of infant-robot interaction profiles. The current
406 GUI and SIPPC 4 robot were tested with 4 infants between the ages of 5 and 9 months, by two therapists
407 over five months at OU Health Sciences Center before deployment to CHOP (the Children's Hospital of
408 Philadelphia) for the clinical study. The test sessions with the infants were conducted using the same clinical
409 study procedures outlined in Prosser et al. [45]. Test sessions consisted of three 5-minute trials of active
410 motor training. Infants were outfitted with the motion capture suit, placed on the SIPPC, and encouraged
411 to move toward toys by the therapists. The therapists tested the data recording, trial control, robot and
412 parameter control, video streaming, sensor calibration, data management, note-taking, robot state display,
413 and infant profile management operations of the GUI. After each testing session, the therapists provided the
414 engineering team with feedback about their experience, including the functionality and ease of use of the
415 GUI and the expected functionality of the SIPPC. This feedback led to adjustments to the organization and
416 operation of the GUI components. In particular, the order of GUI components in the interface was changed
417 to reflect the temporal order in which the therapists tended to interact with these different components. In
418 addition, we adjusted the range of selectable control parameter values to better reflect the useful range of
419 values, as perceived by the therapists.

420 The GUI is also robust to different operating systems (OSs) since web applications are designed to be
421 browser-based and thus largely hardware and OS agnostic. The GUI was tested on Linux, Windows, Android,
422 iPhones, and iPads. The GUI layout and styles dynamically changed to accommodate different screen sizes
423 and resolutions for various devices such as phones and laptops. The GUI was also tested using the Edge,
424 Chrome, Safari, and Firefox web browsers and operated similarly during the sessions in each.

425 **Data Collection Protocol**

426 The SIPPC 4 is one of three stages in a larger, ongoing study on locomotor development of infants at high
427 risk for CP [45]. The first stage is Early Spontaneous Neuromotor Behavior, the second stage is Prone
428 Locomotor Training with the SIPPC, and the third stage is Upright Locomotor Training. The larger study
429 will enroll 60 infants. Presently, over 30 infants have participated in the prone stage. The inclusion criteria

430 for the infants in the larger research study are a brain injury at birth without other known genetic conditions
431 unrelated to CP. Infants are considered to be at high risk for CP if they have Absent Fidgety movements on the
432 General Movements Assessment at 3 or 4 months of age [46] and/or score below the 16th percentile on the Test
433 of Infant Motor Performance (TIMP) Scales [47] at 4 months of age.

434 For the second stage of the study with the SIPPC, infants participate in prone training sessions, with each
435 therapy session ideally consisting of three trials (in practice, the number of trials depends on the disposition
436 of the infant). Active training time during each trial is 5 minutes [45]. Sessions are conducted three times a
437 week starting when the infant is 5 months of age until the infant is 9 months of age (ages are corrected if the
438 infant is born before 37 weeks) or until the infant begins to crawl independently. Prone training sessions occur
439 at the infant's home, at their inpatient bedside, or in a Center for Rehabilitation outpatient location (i.e., CHOP),
440 depending on the most convenient location for the family. Currently data collection is ongoing. Further details
441 on the participants and the study protocol are in [45].

442 For each therapy session, the therapist work with the caregiver to use time with the infant as effectivly
443 as possible. After powering on the SIPPC, a period of four minutes is required to ensure all components
444 are booted and initialized. During this time, the motion capture suit is placed onto the infant before placing
445 and securing the infant prone onto the platform. Additionally, the therapist open the Session View page of the
446 GUI (Figure 6) to verify that all the sensors are connected. Once the SIPPC is ready, the Session View page is
447 used to set up and control the SIPPC to guide motor training. The appropriate infant profile is selected based
448 on the infant code number (prefix "LL0") from a drop-down menu on the Session View page. Before each trial,
449 the suit and force/torque (F/T) sensors are calibrated by pressing the corresponding buttons in the GUI. The
450 initial calibration of these sensors is paramount to the protocol and thus must be completed before trials may be
451 initiated. While the IMUs are strapped to the infant, their orientation relative to the joints is unknown. These
452 orientation offsets are measured by placing the infant into a known configuration, with both arms extended to
453 90 degrees of abduction (to form a 'T'), the trunk horizontal, and the legs straight with the feet plantarflexed.
454 As each segment is positioned, the therapist presses the corresponding calibration button (Figure 6). The F/T
455 sensor is calibrated by pressing the "Force/Torque Calibration" button in the GUI (Figure 6), to tare out the
456 infant's weight and mass distribution.

457 The therapist then checks the three camera views in the GUI's video stream to verify that they capture
458 all the infant's limbs and the objects in front of the infant. The robot response and infant control are engaged
459 before starting the trial timer by using the corresponding buttons in the Trial Control panel of the Session View
460 page (Figure 6). The trial timer will automatically stop once 5 minutes of active prone training has elapsed.
461 The therapist can also pause and resume the trial timer, as needed, to turn the SIPPC away from obstacles
462 or untangle any cabling or hoses that are connected to the infant (these connections are not related to the
463 study). During each trial, the infant is encouraged by the therapist or caregivers to move toward motivating
464 objects (e.g., toys, teethers, or caregivers) and practice lifting their head, shifting weight, reaching, kicking, and
465 other crawling-relevant movements. During these sessions, the therapist uses the GUI (typically the Session
466 View or Parameter Page; Figure 7) for any additional modifications to the robot responsiveness in real-time,
467 based on the infant's performance, skill level, mood, environment, or unique needs for practicing and promoting
468 locomotion.

469 **Results: Therapists-Robot-Infant Interaction**

470 Most data flowing through the SIPPC system are logged for the entire duration of the data collection session;
471 each element of data is time-stamped at a millisecond resolution. Logging commences once the therapist
472 provides key session data, including subject and therapist ID, and continues through all session trials. The
473 logged data include robot state, infant kinematics, all GUI input events, and any response that the system

474 has to these requests from the therapist. This information allows us to assess how the therapist is using the
475 GUI within individual sessions and longitudinally across the full set of sessions. For example, we can measure
476 1) the frequency at which the therapist makes adjustments to the SIPPC control parameters and when these
477 changes are made (specifically, preceding a trial vs during a trial; Figures 8–10), 2) the frequency of minor
478 procedural deviations for situations in which the therapist has some latitude in deciding when or whether to
479 perform certain activities (e.g., whether the suit is re-calibrated between trials; Figure 10), 3) longitudinal
480 changes in the use of the different control modes and associated parameters (Figures 8–11), and 4) the
481 number and degree of completion of each trial within a session.
482

483 Usability

484 The therapists (and sometimes caregivers) were the main users of the GUI. The therapists report that
485 the GUI operated correctly during most sessions. The therapists describe the infants as demonstrating a
486 wide range of motor skills and motivation to move at the start of SIPPC therapy sessions. The therapists
487 report changing the parameters over time to adapt to the infants and their changes when appropriate.
488 Generally, the therapist sets the GRF threshold low in the initial training sessions. This causes more frequent
489 robot movements in response to any movement initiated by the infant, including random weight shifts and
490 spontaneous extremity movements. The discrete movement distance parameter was also set high so that the
491 robot could travel to desired items with few steps.

492 With time, many infants demonstrated more consistent, coordinated movements (e.g., reaching, head
493 lifting, kicking and weight shifting) to move the robot, indicating that the infant was learning how to use
494 their intentional movements to move in the environment. In response, the therapists adjusted the parameters
495 to make it more challenging for the infant by increasing the GRF threshold and decreasing the distance that
496 the robot travels. The infant then had to produce larger forces and move more to reach desired items in the
497 environment.

498 For some infants, initial interactions with the SIPPC, or being prone, required the therapist to use an
499 alternative approach. An example from the therapist's experience is one infant who did not tolerate being
500 prone on the platform and had very limited tolerance to being prone off of the SIPPC. On the SIPPC, this
501 resulted in increased weight shifts, limb movements, and changes in GRFs. The SIPPC's response to these
502 movements and GRF changes was forward motion that further upset the infant. For this infant, being on
503 the SIPPC was unpredictable and resulted in the infant being unable to achieve a calm state for engaging in
504 the training session. Therefore, the treatment goal for this infant began with increasing tolerance for prone
505 positioning on the SIPPC by decreasing movement of the robot. The therapists used the GUI to reduce the
506 robot's response by decreasing the forward distance the robot traveled in response to the infant's movement
507 and to increase the GRF threshold required to move the robot. The gesture recognition modes for the suit
508 were also turned off to reduce the likelihood of unintended robot movements in response to limb movements.
509 This dampened the SIPPC's responses such that if the infant produced spontaneous limb movements, it
510 resulted in minimal response from the SIPPC, creating a more stable and predictable play environment for
511 the infant. As the infant's tolerance for prone positioning improved, the infant was able to be on the SIPPC
512 while remaining calm engaging in play. Then, the therapist was able to adjust the robot parameters to
513 gradually increase the responsiveness of the SIPPC. By the end of the intervention time, the infant was able
514 to tolerate prone positioning for all three trials, but was not able to produce much movement within his
515 environment. Although this infant did not achieve proficiency with moving the SIPPC to interact with the
516 environment by the end of the 16 weeks of training, the system was able to be modified to meet the individual
517 needs and current skill set of the infant and gradually advanced as this skill set improved.

518 In general, the therapists report that the GUI performed consistently across therapy sessions. However,

519 one common critique among the therapists is the inconvenience of scrolling between the video stream and
520 main session operations on the Session View page while working with the infant. One therapist would zoom
521 out on the page to include more components within a single view. Overall, the therapists reported that major
522 issues were rare and found the GUI easy to use and intuitive once they had been trained to use the system.

523 When preparing the SIPPC for a therapy session the therapists report needing to spend 5–10 minutes
524 in both the home and inpatient settings. Calibrating the suit took less than 30 seconds. In addition, the
525 therapists often checked the data streams and sensor status to verify proper functioning of the full system.
526 Latencies in these streams are on the order of 10s to 100 ms and are generally not perceived by the therapist.
527 The GUI is configured to receive robot and infant data at 10Hz, with the video stream averaging 8 to 9FPS
528 during live therapy sessions. In some error conditions (e.g., low data storage space, loose sensor connections),
529 detailed diagnostic information from the GUI was sufficient for the therapist to resolve the problem. In other
530 instances, the therapists would have to restart the system, delaying the start of the session. On rare occasions,
531 some issues such as an interfering WiFi access point or a mechanical failure could not be solved in the field
532 by the therapist.

533

534 **Preliminary Data**

535 In Figures 8, 9, and 10, we illustrate the data that is being captured, present examples of how the therapist
536 uses the GUI, and describe how the therapist changes to control system parameters relate to infant activity.
537 An overview of concurrent infant, robot, and GUI parameter data during a single trial for the infant with the
538 subject code LL011 is shown in Figure 8. A subset of the parameter values and their changes are shown under
539 Therapist Controls in the figure. The times when the infant triggers a discrete assistance and the type (indicated
540 by color), are shown under Robot Assistance. The changes in the robot velocities and detected forces are under
541 Robot State. The speed of the infant’s wrists and ankles and the head height are under Infant Kinematics.

542 In the GUI, therapists can pause trials to temporarily disable the infant’s control of the robot to manually
543 drive the robot and reposition it when the infant runs out of forward-moving space. The first row, labeled as
544 Trial Paused, has red patches indicating the duration and times when the therapists pause the trial (the pause
545 button is found under the Trial Control panel in Figure 6). Under Therapist Controls, in the rows corresponding
546 to Assistive Movement Linear Distance, PS Linear Threshold, and PS Rotational Threshold, we see the values
547 of several control mode parameters that therapists can change to adjust the robot sensitivity. Assistive Linear
548 Distance corresponds to the Distance parameter under Discrete Assistance on the SIPPC Parameters page (Fig-
549 ure 7). Decreasing the distance parameter reduces how far the SIPPC will travel after a discrete assist (e.g.,
550 power steering) is triggered by the infant. The therapist increases the Power steering (PS) Linear Threshold
551 twice, making triggering PS assistance more challenging for the infant. The red bars under Robot Assistance
552 indicate that the infant generated enough GRF to trigger PS assistance. The peaks in GRFs (i.e., Forward
553 Force and Turning Torque) can be observed whenever PS is triggered. The light-blue box starts slightly
554 after 100s. A few seconds before this time frame, the infant triggers left movements using Gesture 1 (G1)
555 assistance, followed by forward movements using PS. These are indicated by the blue bars and the following
556 red bars, respectively, under the Robot Assistance section. Under the Robot State section, there is a corre-
557 sponding positive peak in the rotational robot velocity that occurs just before the left G1 events. Following
558 the left G1 events are the forward PS events, that correspond to a positive square peak in linear robot velocity
559 under the Robot State section. Then, the trial is paused as indicated by the red patches under the Therapist
560 Controls for Trial Paused. After the trial is paused, the therapist manually drives the robot backward and then
561 to the right, as seen by the black bars in the Robot Assistance section. Under the Robot State section, there are
562 small corresponding negative peaks in the linear robot velocity when moving the robot backward and negative
563 peaks in the rotational robot velocity when turning the robot to the right (at around 150s). G1 is triggered by
564 infant limb gestures, therefore, we see peaks in the limb speeds under Infant Kinematics. However, due to the

565 complexity of gestures that trigger robot assistance events, it is difficult from just the limb speeds to describe
566 the precise gesture that caused the assistance event. Another interesting observation from this figure is the head
567 height under the Infant Kinematics section. The head height was relatively constant then started increasing
568 sometime before 100s. The peaks indicate the infant was moving their head up and down. Then before 200s,
569 the head height appears to stabilize.

570 Four example trials for one subject (infant LL013) at different weeks in the study are shown in Figure 9.
571 In particular, it shows longitudinal changes in infant activity, changes in how the infant engages with the
572 SIPPC, and active changes the therapist is making to the control system parameters during the trials. These
573 observations are also supported by Tables 2a and 2b. During weeks 4 and 8 fewer overall Robot Assistance
574 events are triggered than in weeks 12 and 16. The number of PS (red) to G1 (blue) infant-triggered assistance
575 events is about the same in weeks 4 and 8 (Tables 2a and 2b). In weeks 12 and 16, G1 is triggered at least
576 twice as often as PS. The average assistance counts over 4-week periods are shown in Table 2b. A linear
577 model is fit to the 4-week averages in Table 2b and the coefficient of determination (R^2) and its corresponding
578 p-value are computed to identify increases in the use of the robot assistance as the infant matures. There
579 is a significant increase in the average number of G1 events throughout the study ($p < 0.008$). There is no
580 strong trend in the average number of PS events ($p < 0.574$). The average number of PS events increases
581 in the middle of the study and then drops in the last four weeks. Additionally, we observe changes in how
582 the therapist modifies the parameters over time. In week 4, the therapist adjusts parameters multiple times
583 throughout the trial. In week 8, the therapist adjusts the parameters once in the middle of the trial. In weeks 12
584 and 16, the therapist adjusts the parameters within the first 100s of the trial. Lastly, most of the infant-triggered
585 robot assistance events are in the forward direction and the number of events increases from weeks 4 and 8 to
586 weeks 12 and 16.

587 Figure 10, along with Figure 11, illustrate the high-level control mode choices that the therapists are
588 making longitudinally and their potential impact on assistance events across the infants. As part of the therapy
589 standard procedure, therapists calibrate the suit (blue) at the beginning of each trial, though is not required so
590 as to give the therapist flexibility in working with the infant. Out of 1382 trials, the suit is calibrated for 89%
591 of them.

592 Therapists also have the option to write notes (Session Notes; orange) in the GUI during trials. This feature
593 is used on average for 72% of the trials for most infants. Additionally, the therapists used the GUI to enable
594 the control modes Power Steering (PS; green) and G1 (Gesture 1; red) during most trials (99% and 85% of
595 the trials, respectively). In particular, the trials where each control mode is enabled for more than 50% of the
596 trial duration, are seen in Figure 10. G2 (Gesture 2; purple) is typically enabled when G1 is turned off, such
597 as in weeks 4 through 12 for infant LL041. G2 is typically unused during trials early in the study because G2
598 requires more precise movements to be triggered; exceptions are infants LL006 and LL015. Generally, G2
599 is enabled in trials later in the study, such as for infants LL016 and LL028. For infants, LL040-LL042, the
600 therapists enabled G2 for at least half of the weeks in the study. In Figure 11 are the 4-week averages of the
601 assistance counts across infants. As expected, infants who have G2 enabled for a consistent time frame have
602 non-zero average event counts during one or more 4-week periods. G1 appears to have the most counts
603 more often than PS. However, these counts are also impacted by therapist decisions for parameters such
604 as the PS thresholds. Figures 9 and 11, and Table 2 are consistent with the observation for infant LL013
605 that the PS is triggered less than G1. Increases in PS thresholds seen in Figure 9 suggest the influence
606 of these parameters on potentially favoring G1 over PS, since increasing PS thresholds makes triggering
607 PS more challenging. A more detailed analysis comparing kinematics, GRFs, assistance counts, and robot
608 behavior with the therapist's chosen parameter values will provide further insight into the infant behavior and
609 development. It will be discussed more in future work.

Infant LL013	Power Steering	Gesture 1
Week 4	13	11
Week 8	10	15
Week 12	16	32
Week 16	8	32

(a) Discrete assistance counts for trials in Figure 9

Infant LL013	Power Steering	Gesture 1
Week 1 to 4	7.2	15.9
Week 5 to 8	9.1	20.3
Week 9 to 12	9.9	25.7
Week 13 to 16	4.1	33.2
R^2	0.1815	0.9841
p-value	0.574	0.008

(b) 4-week averages of the discrete assistance counts for infant LL013

Table 2: Example of 4-week assistance count data. (a) lists the total number of Power steering (PS) and gesture 1 (G1) events for each trial in Figure 9, and (b) lists the 4-week averages of the number of discrete assistance events for infant LL013. This can provide a general understanding of changes in infant engagement with the discrete assistance events in an environment where the therapist’s parameter decisions also affect the counts. A linear model is fit to the average assistance counts to compute the strength of the relationship between the week in the study and the assistance counts. This is indicated by R^2 (coefficient of determination) and the corresponding p-value. This is an example of one type of analysis that can be performed to understand and summarize longitudinal patterns for each infant.

Discussion

Infants at high risk for CP often have limited opportunities to experience movement than their peers and many do not learn independent prone locomotion skills. The development of the SIPPC 4 assistive crawling robot was guided by prior versions [21–23, 27, 44] and our proposed design principles. The current version of the GUI augments and expands upon the contributions of the prior SIPPC studies. Based on preliminary data and feedback from the therapists, the GUI facilitates personalized motor training opportunities for infants at risk of CP and with motor delays. Prone training on the SIPPC, in conjunction with guidance from the therapists using the GUI, provides infants with CP with more opportunities to practice prone skills that will encourage locomotion (e.g., head lifting, weight bearing on the arms, and reaching) while teaching the infants that purposefully initiated actions can result in successful movement towards a desired goal. Using the GUI to adjust the feedback also offers an intuitive method of communication with the learner, particularly when communication challenges are present (e.g., with infants).

The GUI provides an approach for data monitoring and documentation to observe therapist-robot-infant interaction and infant performance (Figures 8–12). The unique design of the SIPPC GUI as a web application makes it uniquely available to all authorized users with any WiFi-enabled devices (e.g., laptop, tablet, phone) in their favorite browser. No special installations or configurations are required for therapists or caregivers to easily access the GUI to monitor infant and robot behavior or modify the therapy. The GUI design also makes it easy for authorized users to share data and parameter configurations with post-processing systems or between other SIPPC robots, all from the same interface.

From Figure 9 and Tables 2a and 2b, Infant LL013 showed an increase in Gesture 1 (G1) movements, suggesting potential improvements in limb movement generation and coordination, as well as influences from the therapist’s parameter decisions. The Power steering (PS) thresholds may be high enough such that G1 assistance events are more likely to trigger instead (since only one discrete movement may be triggered at one time). Data available through the GUI provides the potential to monitor longitudinal changes in locomotor behavior, skill, and their relationship to actions taken by therapists to guide motor training. A variety of control modes (Session Configuration section in Figures 6 and 10), with variable levels of assistance

636 (Parameters in Figures 7 and 8), are available for adjustment by the therapist in the GUI. The control modes and
637 assistance levels allow infants at any developmental stage and motor level to explore and practice actions that
638 are precursors to crawling (design principle 1). For example, the infant who the therapists reported initially did
639 not tolerate being prone. The therapists were able to change the parameters in the GUI and modify the robot
640 response such that the infant could acclimate to being prone and on the SIPPC. From Figure 10, although we
641 observe the therapists enable PS more often than G1 for most sessions (99% and 85% of trials, respectively),
642 we often see more G1 assistance events in many infants. This further implicates the need to analyze therapist
643 treatment decisions with infant performance and robot response. This will allow us to quantitatively describe
644 infant-specific choices therapists make. Such analyses will be performed in future work after completion of
645 data collection.

646 From Figure 10, there are some consistent aspects of the therapy implementation, such as an adherence to
647 the suit calibration before most trials. For approximately 89% of trials, the suit is calibrated before beginning
648 subsequent trials. The GUI only mandates the suit calibration for the first trial at the start of the session
649 and provides flexibility in calibrating subsequent trials for therapists to accommodate unique infant needs. In
650 addition to suit calibration, the therapists choose to enable PS and G1 most of the time and consistently
651 disable G2 early in the study for most infants. If G2 provides a useful training strategy early on for an infant,
652 such as LL041 and LL042, then the therapists continue using G2 throughout most or all of the remainder of
653 the study. Since the robots can share subject parameter configurations, this can easily be communicated to
654 all therapists conducting therapy sessions for a particular infant.

655 Motivated by design principle two, the robot acts as a mediator from performance to results. The robot
656 measures the infant's real-time performance using the kinematics and GRFs, and augments the infant's actions
657 to encourage continued motor exploration. The infant performance relative to the robot response can be seen
658 in Figures 8 and 9. We can identify when the robot moves and which limbs contribute to the movement.
659 For example, in Figure 8, around 0s, a left robot movement is triggered due to G1 (blue bars under Robot
660 Assistance). The peaks in the Right Wrist Speed and the Left Foot Speed suggest they are the likely contributors
661 to triggering the assistance from G1. There is also a peak in the Robot Rotational Velocity, which is a left turn
662 of the robot. At around 100s, peaks in the forward force trigger power steering and a forward robot movement.
663 Most of the limbs likely contributed to the generation of this net force since multiple and frequent changes in
664 limb speeds are observed under the Infant Kinematics. Additionally, in this particular trial, the right foot sensor
665 was not collecting data due to a damaged cable. However, the SIPPC and GUI still operate functionally for
666 therapeutic purposes since the SIPPC has control modes that are independent of the suit (i.e., Force Control
667 and Power Steering). This suggests that multiple control strategies defined by different performance measures
668 (e.g., force vs gesture) can be an approach to help enhance robustness to failure in clinical settings and reduce
669 the likelihood of missed or failed sessions.

670 In Figure 9, we observe potential longitudinal behavioral changes in an infant LL013 as the changes in
671 forces and limb speeds. We also see changes in the relative adjustments made by the therapists to the robot
672 parameters (e.g., changes in PS thresholds). From the first 4 weeks to the last 4 weeks, the infant is triggering
673 the robot assistance more frequently on average (Table 2b). G1 is triggered relatively more than PS in the
674 later weeks. This might be due to changes in the limb movement and coordination relative to changes in
675 strength to generate GRFs. However, this will also be impacted by the PS Thresholds. For example, the PS
676 Linear Threshold is increased by the therapist near the beginning of the trials in weeks 12 and 16, whereas, the
677 threshold is increased later in the trials in weeks 4 and 8. This increase in the PS threshold makes it more
678 challenging to trigger PS and potentially increases the likelihood for G1 to be triggered instead since only one
679 discrete assistance can be triggered at a time. This behavior could be ideal when an infant has met sufficient
680 strength goals and the therapist wants to encourage the generation of more coordinated limb movements.

681 Moreover, the real-time robot response adjustment capabilities provided through the GUI make it possible
682 for the therapists to help address real-time changes in infant mood and behavior as well as longitudinal changes

683 in motor skills. For instance, from the therapist reports in the Usability Section, therapists can support infants
684 who are initially responsive to therapy and infants who initially don't respond well to various aspects of prone
685 training. The therapists were able to provide and progressively increase the challenge of the therapy session
686 for the infants that were open to being on the SIPPC. As these infants became advanced in generating
687 goal-directed robot movements (e.g., continuous movement towards toys), the therapist can immediately
688 reduce the sensitivity of the robot response by either removing the suit recognition modes, increasing the
689 GRF threshold for PS, or reducing the distance that the robot travels in response to infant movements. In
690 situations where infants are uncomfortable on the SIPPC, the therapists also reported being able to guide
691 these infants toward obtainable motor and behavioral goals. This was typically done by reducing the sensitivity
692 of the robot to avoid upsetting the infant. This was done by disabling the gesture recognition modes and
693 reducing the distance the robot moves. Once the infant became comfortable on the SIPPC, these control
694 parameters were gradually and slowly restored to typical starting values.

695 Consistent with design principle three, the SIPPC GUI provides the therapist with detailed feedback (Fig-
696 ures 6, 8–12) that can be used dynamically to modify the behavior of the robot to further guide locomotor
697 learning. By adjusting the robot behavior, motor training difficulty can either be decreased to motivate the in-
698 fant or increased to challenge the infant. The ability to adjust the difficulty is crucial to prevent dependency on
699 the robot-assistance, ensure that progress translates to daily life [4, 16] and is necessary to provide personalized
700 care for precision rehabilitation [19]. This approach can provide the infants with more opportunities to move
701 within their abilities and skill level while also allowing them to practice increasingly difficult actions.

702 **Limitations**

703 While many useful outcomes have been observed with the current application of the SIPPC-GUI system, some
704 challenges have also been encountered. The current analyses are limited quantitatively because data collection
705 is still ongoing. In future work, the complete analysis will compare all data from the SIPPC-GUI system (i.e.,
706 robot velocities, robot forces, robot parameters, infant kinematics, force, and video data, and GUI events).
707 Presently, we describe the GUI design, its behavior and function, and a detailed qualitative analysis of the
708 preliminary data to lay the groundwork for a more comprehensive analysis that will be done as future work.

709 Another challenge is the 4–5 minutes start up time for the robot and GUI to finish booting and initializing all
710 components. This is particularly problematic as some infants have limited tolerance for waiting for a therapy
711 session to start. Because of the interplay between the large number of hardware and software components,
712 there exist substantial constraints on this boot-up sequence. Current timing was selected conservatively to
713 ensure that all these timing dependencies were met. This full boot-up process will be evaluated in future
714 work as the current version of the robot is mainly for research use.

715 Additionally, when therapists pause the infant's control to manually drive the robot, there are some instances
716 when the infant is still performing an action. An example of this can be seen in Figure 8 after 100s. The infant
717 is still moving the limbs at similar speeds while the therapist has paused the trial. Due to the trial pause, this
718 action will have no impact on the robot response, even for intentional movements. Therefore, these pauses
719 should be as short as possible, especially when the infant becomes more skilled.

720 Another limitation of the system is related to the camera view alignment. Camera alignment is not automatic
721 and must be adjusted manually by the therapist such that the infant's limbs are clearly visible in the video feed
722 shown in the GUI. Failing to take this step can negatively impact the human-driven analysis of the infant-robot
723 behavior. Future system components could automatically analyze the incoming images and confirm that the
724 alignment is sufficient.

725 The system can also be impacted by certain physical attributes of the therapy environment. In particular,
726 certain flooring, such as some tiles and carpets, can alter the mobility of the robot, whereas infants might have
727 little issue crawling over these same flooring conditions. Therefore, the parameter settings might need to be

728 adjusted more often early in a session to improve the robot response to the flooring. The last limitation of the
729 system is that the GUI presents a fixed organization of the graphical components. Future work could allow the
730 therapists to customize the layout of Session View to make their interaction with the GUI more efficient and
731 specific to their particular style of interaction.

732 **Clinical Implications**

733 The SIPPC offers a flexible and transportable system to adapt to the unique developmental path of each infant.
734 It is crucial to provide rehabilitation that meets the individual needs of each infant. The design principles we
735 have proposed in conjunction with the SIPPC and its GUI enable therapists to customize therapy contingent
736 upon the specific needs of the learner and adjust the motor training difficulty to the unique skill level of each
737 individual. In particular, infants with CP can face various combinations of challenges such as reduced strength,
738 endurance, or coordination. The selection and configuration of GRF-based control modes such as PS, can
739 address strength and endurance training. The gesture-based control modes (G1 and G2) can support infants
740 who struggle with coordination or strength. Moreover, the combination of these types of control modes can
741 target all these challenges for an infant who needs practice with each of them. The GUI is a valuable tool
742 in providing a unique and individually tailored training program while utilizing the same general intervention
743 approach for a variety of infant skill levels.

744 When providing intervention to facilitate developmental skill acquisition, it is important to provide a train-
745 ing environment that includes skill specificity, motivation for the infant, and modulated activity appropriate to
746 the current skill level [7, 13, 21]. Using the GUI with SIPPC training allows for all of these conditions to be
747 met and applies to a variety of infants. The SIPPC allows for crawling-like movements to result in forward
748 motion toward desired toys (Figure 12). The GUI allows the activity to be easily and efficiently upgraded or
749 downgraded to meet an appropriate training intensity for each infant (Figures 6–10). This feature is particularly
750 important to tap into motor learning principles for developmental skill progression and precision rehabilitation.
751 The activity needs to be challenging enough to continue to entice the infant into performing new and increas-
752 ingly difficult skills, without being too challenging and resulting in decreased motivation for task performance.

753 **Future Work**

754 Future work will aim to improve the therapist’s interaction with the robot through the GUI. Improving the boot-
755 up speed of the SIPPC and GUI can further optimize time spent in motor training with infants. Additionally,
756 wider angle video feeds or automated positional control of the cameras in the GUI can further increase the time
757 that infants spend in motor training with the guidance of the therapist. However, the cameras are predominately
758 for research use and not for pure clinical applications. The large amount of data available in the GUI can be
759 overwhelming, therefore the ability of the therapists to quickly and easily customize the data presented in the
760 GUI can further improve the quality of prone locomotor training for infants.

761 Future quantitative longitudinal analyses of the concurrent infant, robot, and parameter session data will be
762 used to measure the design performance. For example, coders evaluate infant behavior and motor skills from
763 the video collected during the sessions. These data will be compared with robot and parameter data to assess
764 how therapist control changes impact infant locomotor behavior.

765 Future work should also apply these design principles to other assistive therapy devices, age groups, and
766 motor disorders. The general design and layout of the GUI, and overall SIPPC system, can be adapted to other
767 groups and devices such as robotic assistive gait devices. In the current system, user actions in the GUI send
768 change requests to the robot control system. This same workflow can be applied to the corresponding control
769 system of any robotic device and provide therapists with detailed control over the behavior of the robotic
770 therapy device. Additionally, ROS is used to manage data communications between the robot controls, the

771 sensors, and the GUI. This can be applied to the integrated sensors that are a part of any robotic device
772 to send sensor data to the GUI for display to the user. These data can be displayed in the manner most
773 pertinent to the therapy objectives.

774 **Acknowledgements**

775 This work was funded by the NIH/NICHD IR01HD09836Y Grant. Research reported in this publication was
776 also supported by the Eunice Kennedy Shriver National Institute Of Child Health & Human Development of
777 the National Institutes of Health under Award Number R01HD098364. The content is solely the responsibility
778 of the authors and does not necessarily represent the official views of the National Institutes of Health. We
779 would like to thank Caroline A. Wagner for her work during test sessions with the infants, SIPPC and the GUI,
780 Reza Torbati for his contributions to the video compositing software, Hope Chambers for her contributions to
781 fabrication efforts, Zachary Muller for his contributions to the software and data analysis methods, and Farzana
782 Tasnim for her contributions to the data analysis methods.

783 **Conflicts of Interest**

784 None

785




Figure 4: Data Flow Diagram. Active communications, in the rounded rectangles, refer to data transmissions during therapy sessions, such as the movement of video data between the main and video RPi's. Inactive communications, in the rectangles, involve moving data before or after sessions. Recorded session and robot data stored on the robot are uploaded to secure remote cloud storage and session configurations, parameters and notes can be downloaded, allowing different robots to be interchangeable. The control state is comprised of sensor data from the force/torque (f/t) sensor and information about the wheels' state. The control state also provides the connection status of the f/t sensor and the wheels. Control state change requests are petitions to alter the state of the wheels. The suit state consists of sensor data from the IMUs in the motion capture suit, in addition to the connection status of each sensor. The camera state contains the parameters brightness, contrast, exposure, and focus values for all three cameras. A change request for the camera petitions to adjust any of the camera parameters. The video stream is a composited video of the three camera views and a set of data plots. The session state is an amalgamation of information about the current session and trial, such as the subject number, if data recording has started, if the trial has started, the trial number, and the attending therapist. Session state requests petition to alter aspects of the session state such as requesting to start recording or selecting a subject. Processors states describe the CPU load and available memory and disk space for each of the RPi's. The recorded robot, video, and session data are uploaded to secure cloud storage after the session has ended and once the robot is connected to the internet via an Ethernet connection on the main RPi. The subject parameter sets and session notes can be uploaded and downloaded as well.





Figure 5: Diagram of network and device connections. The robot has five specialized processors (RPi's and two Mbed microprocessors) that communicate via Ethernet. The main RPi acts as an access point and web server broadcasting the SIPPC local network and GUI. The robot is connected to the internet through main RPi to synchronize data between secure remote machines and other SIPPC robots. The video RPi (sippcvideo) is dedicated to capturing and compositing the camera videos and data plots. The control Mbed (sippccontrol) monitors and manipulates the wheels and force/torque (F/T) sensor. The suit Mbed is connected to 12 IMUs on the motion capture suit.




Figure 6: Session view is the main page, consisting of the most relevant operations necessary for the therapists to perform therapy sessions and provide the infant with as much attention as possible. Therapists receive a summary of the robot status (System Status), can select the infant code and therapist name, calibrate the suit, view velocity, force, and infant limb data in real-time, modify robot control modes and parameters (Session Configuration), manage the trial, drive the robot (Manual Assist), and write notes about the active session.

SIPPC Parameters

[Save Configuration](#) [Export Parameters](#)

Discrete Assistance

Distance (cm)	Turning (deg)
<input type="text" value="20"/>	<input type="text" value="45"/>
< 0 <input type="range" value="20"/> 2 >	< 0 <input type="range" value="45"/> 2 >
Duration (sec)	Refractory Time (sec)
<input type="text" value="1.25"/>	<input type="text" value="0.2"/>
< 0 <input type="range" value="1.25"/> 2 >	< 0 <input type="range" value="0.2"/> 0.5 >

Power Steering

Forward Thres (N)	Turning Thres (Nm)
<input type="text" value="0.5"/>	<input type="text" value="0.4"/>
< 0 <input type="range" value="0.5"/> 4 >	< 0 <input type="range" value="0.4"/> 2 >

Force/Torque Sensor

τ (sec)	Torque Filter Scale
<input type="text" value="0.2"/>	<input type="text" value="0.5"/>
< 0.01 <input type="range" value="0.2"/> 2 >	< 0.01 <input type="range" value="0.5"/> 1.5 >
τ Tare (sec)	Auto Tare
<input type="text" value="60"/>	<input type="checkbox"/>
< 1 <input type="range" value="60"/> 100 >	

Force/Torque Control

Global Force/Torque Gain	Backward Force Boost (N)
<input type="text" value="1"/>	<input type="text" value="0"/>
< 0 <input type="range" value="1"/> 2 >	< 0 <input type="range" value="0"/> 4 >
Forward Force Boost (N)	Forward Force Damping (N)
<input type="text" value="1"/>	<input type="text" value="1.5"/>
< 0 <input type="range" value="1"/> 4 >	< 0 <input type="range" value="1.5"/> 5 >
Forward Force Deadband (N)	Forward Force Max (N)
<input type="text" value="1"/>	<input type="text" value="0.2"/>
< 0 <input type="range" value="1"/> 3 >	< 0.2 <input type="range" value="0.2"/> 4 >
Turning Torque Boost (Nm)	Turning Torque Damping (Nm)
<input type="text" value="2"/>	<input type="text" value="3"/>
< 0 <input type="range" value="2"/> 4 >	< 0 <input type="range" value="3"/> 10 >
Turning Torque Deadband (Nm)	Turning Torque Max (Nm)
<input type="text" value="0.02"/>	<input type="text" value="0.5"/>
< 0 <input type="range" value="0.02"/> 0.3 >	< 0.2 <input type="range" value="0.5"/> 2 >

Figure 7: Parameter page for adjusting how the robot responds to infant movements. Discrete assistance parameters specify the distance, angle, duration, and time between discrete robot assistance triggers. Power steering parameters specify the threshold GRFs that trigger robot assistance. Force/Torque Sensor parameters manage the filtering of the sensor noise. Force/Torque Control parameters specify the magnitude of the proportional force assistance.







Figure 8: Summary of select robot and parameter data during one trial for infant LL011, during Week 4. The section titled Therapist Controls contains some of the parameters in the GUI that the therapist can adjust to personalize the motor training trial (PS is power steering). Therapists can pause the trial (Trial Paused), adjust the linear distance the robot moves when the infant triggers an assistive response (Assistive Movement Linear Distance in centimeters), adjust the linear GRF threshold (in Netwons) to trigger PS, and adjust the rotational GRF threshold (in Newton-meters) to trigger PS. The Robot Assistance section shows the discrete robot assistance events that are triggered by the infant’s movements as vertical tick marks. Black indicates manual assistance where the therapist drives the robot, red indicates power steering, and blue indicates Gesture 1 (G1). The Robot State section shows the robot velocities (in cm/s) and sensed forward forces (in Newtons) and turning torque (in Newton-meters). The Infant Kinematics section shows the infant limb speeds (in cm/s) for the wrists and feet, as well as the height of the head (in centimeters). The light blue box highlights a region when the therapist paused the trial, right after the infant made a series of forward movements triggered by Power steering in red under Robot Assistance. While the trial is paused, the therapist manually moved the SIPPC, indicated in black under the Robot Assistance. In this particular trial, the right foot sensor was not collecting data due to wear-and-tare on the cable. However, the SIPPC and GUI still operate functionally for therapeutic purposes despite a sensor missing.




(a) Week 4



(b) Week 8



(c) Week 12



(d) Week 16

Figure 9: Longitudinal changes observed during trials for one infant

GUI and Control Mode Use




Figure 10: Use of GUI operations (suit calibration and session notes) and control modes (Power Steering, Gesture 1, and Gesture 2) across trials and infants. The bars for suit calibration indicate the suit was calibrated before starting the trial. The bars for the session notes indicate that the therapists wrote notes during the session. The bars for the control modes indicate when the therapists enabled the control mode for at least 50% of the trial duration.




Figure 11: 4-week averages of the assistance counts for 19 infants. Assistance counts for Power steering (PS), Gesture 1 (G1), and Gesture 2 (G2) are presented here as averages of every 4 weeks. G2 is rarely turned on, thus few assistance events are triggered for it overall. G1 is more commonly the most triggered assistance type for most infants, such as LL011 and LL013. For some infants, we see a steady increase in the average assistance counts for PS or G1. For some infants, these trends appear relatively constant or decrease. These counts are influenced both by infant behavior and therapist decisions.




Figure 12: Real-time Video Feed. The top two panels capture all the limbs from the tip of the toes to the hands. The bottom left panel captures the front of the infant and what the infant is looking at. The bottom right quadrant shows a subset of the captured data. The kinematic skeleton of the limb and trunk positions are displayed along with the corresponding force/torque, robot velocity, and discrete assistance.

SIPPCC9

Session View Parameters Post Processing Subjects Clinicians User Guide System Details Admin Dashboard

Post Processing Procedures

1. Upload and Sync Config, Video, and Log Files

(a) Subject Configuration

Files 0 / 0

Start

(b) Compress Videos

Files 0 / 0

Start

(c) Videos and Logs

Files 0 / 11

In Progress

2. Delete Old Files

(a) Raw Video

Files 0 / 0

Start

(b) Remaining Files

Files: 7 old / 384 total

Start

Figure S1: Post Processing page for synchronizing data with the cloud database and compressing the video files before uploading. The green WiFi symbol at the top right of the page indicates whether there is an internet connection to perform the data uploading. Data can be deleted only after it has been uploaded to free up space on the robot. The data management steps are presented in the order of their intended operation and cannot be performed while the data is recording for a session.

SIPPC9 Session View Parameters Post Processing Subjects Clinicians User Guide System Details Admin Dashboard

System Status

✔ Communications
✔ Force Sensors
✔ Wheel Left
✔ Wheel Right
✘ Suit

Suit Details

Limb Calibration

✘ Trunk
✘ Left Arm
✘ Left Leg
✘ Right Arm
✘ Right Leg

Suit Sensor Validation

✘ Base
✘ Back
✘ Left Upper Arm
✘ Left Fore Arm
✘ Right Upper Arm
✘ Right Fore Arm

✘ Left Thigh
✘ Left Calf
✘ Left Foot
✘ Right Thigh
✘ Right Calf
✘ Right Foot

pi0 - sippcbridge

- Temperature 55°C
- Load 5.90%
- Disk 91.28 GB
- FreeMemory 6.49 GB

pi1 - sippcvideo

- Temperature 47°C
- Load 4.90%
- Disk 102.26 GB
- FreeMemory 375.00 MB
- DiskData 744.30 GB

pi2 - sippcaux

- Temperature 54°C
- Load 23.80%
- Disk 102.19 GB
- FreeMemory 7.16 GB

SHUTDOWN

Figure S2: System Details page for diagnosing issues with various aspects of the robot system. Users can view the status of each wheel and specific sensors. Users can also view details about the three main robot processors (sippcbridge, sippcvideo, and sippcaux).

References

- 788 [1] M Bax, M Goldstein, P Rosenbaum, A Leviton, N Paneth, B Dan, B Jacobsson, D Damiano, and Ex-
789 ecutive Committee for the Definition of Cerebral Palsy. Proposed Definition and Classification of
790 Cerebral Palsy. *Developmental Medicine and Child Neurology*, 47(8):571–576, 2005. doi: <https://doi.org/10.1017/s001216220500112x>.
791
- 792 [2] Novak I., Morgan C., Adde L., Blackman J., Boyd R. N., Brunstrom-Hernandez J., Cioni G., Damiano
793 D., Darrah J., Eliasson A. C., de Vries L. S., Einspieler C., Fahey M., Fehlings D., Ferriero D. M., Fetters
794 L., Fiori S., Forsberg H., Gordon A. M., Greaves S., Guzzetta A., Hadders-Algra M., Harbourne R.,
795 Kakooza-Mwesige A., Karlsson P., Krumlinde-Sundholm L., Latal B., Loughran-Fowlds A., Maitre N.,
796 McIntyre S., Noritz G., Pennington L., Romeo D. M., Shepherd R., Spittle A. J., Thornton M., Valentine
797 J., Walker K., White R., and Badawi N. Early, accurate diagnosis and early intervention in cerebral palsy:
798 Advances in diagnosis and treatment. *JAMA Pediatrics*, 171, 2017. doi: 10.1001/jamapediatrics.2017.
799 1689.
- 800 [3] D. R. Patel, M. Neelakantan, K. Pandher, and J. Merrick. Cerebral palsy in children: a clinical overview.
801 *Translational pediatrics*, 9(1):S125–S135, 2020. doi: <https://doi.org/10.21037/tp.2020.01.01>.
- 802 [4] Diane L Damiano. Activity, Activity, Activity: Rethinking Our Physical Therapy Approach to Cerebral
803 Palsy. *Physical Therapy*, 86(11):1534–1540, 2006. doi: <https://doi.org/10.2522/ptj.20050397>.
- 804 [5] Catherine Morgan, Johanna Darrah, Andrew M Gordon, Regina Harbourne, Alicia Spittle, Robert John-
805 son, and Linda Fetters. Effectiveness of Motor Interventions in Infants with Cerebral Palsy: A Sys-
806 tematic Review. *Developmental Medicine and Child Neurology*, 58(9):900–909, 2016. doi: <https://doi-org.ezproxy.lib.ou.edu/10.1111/dmcn.13105>.
807
- 808 [6] Diane L. Damiano. Rehabilitative Therapies in Cerebral Palsy: the good, the not as good, and the possible.
809 *Journal of Child Neurology*, 24:1200–1204, 2009. doi: <https://doi.org/10.1177/0883073809337919>.
- 810 [7] C. Morgan, I. Novak, R. C. Dale, Andrea Guzzetta, and Nadia Badawi. Game (goals - activity - motor
811 enrichment): protocol of a single blind randomised controlled trial of motor training, parent education and
812 environmental enrichment for infants at high risk of cerebral palsy. *BMC Neurology*, 14(1):203, 2014.
813 doi: <https://doi.org/10.1186/s12883-014-0203-2>.
- 814 [8] Karen E. Adolph. Learning to Move. *Current Directions in Psychological Science*, 17(3):213–218, 2008.
815 doi: <https://doi.org/10.1111/j.1467-8721.2008.00577.x>.
- 816 [9] Max Horovitz and Johnny L. Matson. Developmental Milestones in Toddlers with Atypical Develop-
817 ment. *Research in Developmental Disabilities*, 32(6):2278–2282, 2011. ISSN 0891-4222. doi: <https://doi.org/10.1016/j.ridd.2011.07.039>. URL <https://www.sciencedirect.com/science/article/pii/S0891422211002964>.
818
819
- 820 [10] D.I. Anderson, J.J. Campos, D.C. Witherington, A. Dahl, M. Rivera, M He, I Uchiyama, and M Barbu-
821 Roth. The Role of Locomotion in Psychological Development. *Frontiers in Psychology*, 4:440, 2013.
822 doi: <https://doi.org/10.3389/fpsyg.2013.00440>.
- 823 [11] Karen E Adolph, Beatrix Vereijken, and Mark A Denny. Learning to Crawl. *Child Development*, 69:
824 1299–1312, 1998. doi: <https://doi.org/10.1111/j.1467-8624.1998.tb06213.x>.

- 825 [12] Q. L. Xiong, X. Y. Wu, Y. Liu, C. X. Zhang, and W. S. Hou. Measurement and analysis of human
826 infant crawling for rehabilitation: A narrative review. *Frontiers in Neurology*, 12, 2021. doi: <https://doi.org/10.3389/fneur.2021.731374>.
827
- 828 [13] A. Gonzalez, L. Garcia, J. Kilby, and Peter McNair. Robotic devices for paediatric rehabilitation: a
829 review of design features. *BioMedical Engineering OnLine*, 20(89), 2021. doi: <https://doi.org/10.1186/s12938-021-00920-5>.
830
- 831 [14] R. Kabir, M. S. H. Sunny, H. U. Ahmed, and M. H. Rahman. Hand Rehabilitation Devices: A Compre-
832 hensive Systematic Review. *Micromachines. Current physical medicine and rehabilitation reports*, 13
833 (7):1033, 2022. doi: <https://doi.org/10.3390/mi13071033>.
- 834 [15] Alan W Salmoni, Richard A Schmidt, and Charles B Walter. Knowledge of Results and Motor Learning:
835 A Review and Critical Reappraisal. *Psychological Bulletin*, 95(3):355–386, 1984. doi: <https://doi.org/10.1037/0033-2909.95.3.355>.
836
- 837 [16] Tomoko Kitago and John W Krakauer. Motor Learning Principles for Neurorehabilitation. *Handbook of*
838 *Clinical Neurology*, 110:93–103, 2013. doi: <https://doi.org/10.1016/b978-0-444-52901-5.00008-3>.
- 839 [17] Danielle Levac, Laurie Wishart, Cheryl Missiuna, and Virginia Wright. The Application of Motor Learn-
840 ing Strategies within Functionally Based Interventions for Children with Neuromotor Conditions. *Pedi-*
841 *atric Physical Therapy*, 21(4):345–355, 2009. doi: <https://doi.org/10.1097/PEP.0b013e3181beb09d>.
- 842 [18] Valerie Pomeroy, Salvatore M Aglioti, Victor W Mark, Dennis McFarland, Cathy Stinear, Steven L Wolf,
843 Maurizio Corbetta, and Susan M Fitzpatrick. Neurological Principles and Rehabilitation of Action Dis-
844 orders: Rehabilitation Interventions. *Neurorehabilitation and Neural Repair*, 25(5):33S–43S, 2011. doi:
845 <https://doi.org/10.1177/1545968311410942>. URL [https://www.ncbi.nlm.nih.gov/pmc/articles/](https://www.ncbi.nlm.nih.gov/pmc/articles/PMC4139494/)
846 [PMC4139494/](https://www.ncbi.nlm.nih.gov/pmc/articles/PMC4139494/).
- 847 [19] Catherine P. Adans-Dester, Nicolas Hankov, Anne T. O’Brien, Gloria P. Vergara-Diaz, Randie M. Black-
848 Schaffer, Ross D. Zafonte, Jennifer Dy, Sunghoon Ivan Lee, and Paolo Bonato. Enabling Precision Re-
849 habilitation Interventions using Wearable Sensors and Machine Learning to Track Motor Recovery. *NPJ*
850 *Digital Medicine*, 3 1:121, 2020. doi: <https://doi.org/10.1038/s41746-020-00328-w>.
- 851 [20] A Moinuddin, A Goel, and Y Sethi. The Role of Augmented Feedback on Motor Learning: A Systematic
852 Review. *Cureus*, 13(11), 2021. doi: <https://doi.org/10.7759/cureus.19695>.
- 853 [21] Thubi H. A. Kolobe and Andrew H. Fagg. Robot Reinforcement and Error-Based Movement Learning in
854 Infants With and Without Cerebral Palsy. *Physical Therapy*, 99(6):677–688, 2019. doi: <https://doi.org/10.1093/ptj/pzz043>.
855
- 856 [22] H.A. Kolobe, A. Fagg, Peter Pidcoe, and P. Williams. Effectiveness of Reward-and Error-Based Move-
857 ment Learning in Enhancing Self-Initiated Prone Locomotion in Infants with or at Risk for Cerebral Palsy.
858 *Physiotherapy*, 101:e781–e782, 05 2015. doi: <https://doi.org/10.1016/j.physio.2015.03.3658>.
- 859 [23] Mustafa A Ghazi, Michael D Nash, Andrew H Fagg, Lei Ding, Thubi H A Kolobe, and David P Miller.
860 *Novel Assistive Device for Teaching Crawling Skills to Infants*, volume 113, pages 593–605. Springer
861 International Publishing, 2016. doi: https://doi.org/10.1007/978-3-319-27702-8_39.
- 862 [24] Tummy Time Toys | Best Crawling Toys | Crawligator. <https://www.thecrawligator.com/>. Accessed:
863 March 11, 2024.

- 864 [25] Jack Mcnamara. Redbarn Enterprises Inc. <https://redbarn-enter.com/>. Accessed March 11, 2024.
- 865 [26] J. Laut, M. Porfiri, and P. Raghavan. The Present and Future of Robotic Technology in Rehabilitation. *Current Physical Medicine and Rehabilitation Reports*, 4(4):312–319, 2016. doi: <https://doi.org/10.1007/s40141-016-0139-0>.
- 866
- 867
- 868 [27] David Miller, Andrew Fagg, Lei Ding, Thubi Kolobe, and Mustafa Ghazi. Robotic Crawling Assistance for Infants with Cerebral Palsy. 01 2015.
- 869
- 870 [28] Johanna Renny Octavia, Jorn Ockerman, Jamil Joundi, Arno Penders, and Jelle Saldien Lynn Bar-On. Matti: Tangible user interface for engaging patients in physical therapy towards a motivating rehabilitation. *International Journal of Technology*, 15(3):697–708, May 2024. ISSN 2087-2100. doi: <https://doi.org/10.14716/ijtech.v15i3.6118>.
- 871
- 872
- 873
- 874 [29] Luke A Tucker, Ji Chen, Lauren Hammel, Diane L Damiano, and Thomas C Bulea. An open source graphical user interface for wireless communication and operation of wearable robotic technology. *Journal of Rehabilitation and Assistive Technologies Engineering*, 7:2055668320964056, 2020. doi: 10.1177/2055668320964056. URL <https://doi.org/10.1177/2055668320964056>. PMID: 33403122.
- 875
- 876
- 877
- 878 [30] R. Young, K. Sage, D. Broom, Andrew Hext, Nicky Snowdon, and Christine Smith. Evaluating the usability of a co-designed power assisted exercise graphical user interface for people with stroke. *Journal of NeuroEngineering and Rehabilitation*, 20(95), 2023. doi: <https://doi.org/10.1186/s12984-023-01207-7>.
- 879
- 880
- 881 [31] Alexander Wilkinson, Michael Gonzales, Patrick Hoey, David Kontak, Dian Wang, Noah Torname, Amelia Sinclair, Zhao Han, Jordan Allspaw, Robert Platt, and Holly Yanco. Design guidelines for human–robot interaction with assistive robot manipulation systems. *Paladyn*, 12(1):392–401, 2021. doi: [doi:10.1515/pjbr-2021-0023](https://doi.org/10.1515/pjbr-2021-0023). URL <https://doi.org/10.1515/pjbr-2021-0023>.
- 882
- 883
- 884
- 885 [32] J. B. Southerland. Activity recognition and crawling assistance using multiple inexpensive inertial measurement units. Master’s thesis, School of Computer Science, University of Oklahoma, 2012.
- 886
- 887 [33] Daniel Roetenberg, Henk Luinge, and Per Slycke. sens mvn: Full 6dof human motion tracking using miniature inertial sensors. Technical report, XSENS Technologies, April 2013.
- 888
- 889 [34] Mustafa A. Ghazi, Monique O. Shotande, Thubi H. A. Kolobe, Reza J. Torbati, Julie Skorup, Stephanie Walters, Laura A. Prosser Morgan Alcott, Beth A. Smith, and Andrew H. Fagg. Robotic rehabilitation device for promoting prone locomotion skills in infants with cerebral palsy. *IEEE Transactions on Neural Systems and Rehabilitation Engineering*, in preparation.
- 890
- 891
- 892
- 893 [35] ROS. Robot operating system (ros), 2007. URL <https://www.ros.org/>. Last accessed 1 July 2022.
- 894
- 895 [36] Christopher Crick, Graylin Jay, Sarah Osentoski, Benjamin Pitzer, and Odest Chadwicke Jenkins. *Ros-bridge: ROS for Non-ROS Users*, volume 100, pages 493–504. Springer International Publishing, 2017. doi: https://doi.org/10.1007/978-3-319-29363-9_28.
- 896
- 897 [37] Douglas Crockford. Json, 2022. URL <https://www.json.org/json-en.html>. Last accessed 10 July 2022.
- 898
- 899 [38] Raspberry Pi Foundation. Raspberry pi, 2012. URL <https://www.raspberrypi.org/>. Last accessed 10 July 2022.
- 900

- 901 [39] Pallets Team. Flask, 2010. URL <https://flask.palletsprojects.com/en/2.1.x/>. Last accessed
902 10 July 2022.
- 903 [40] Benoit Chesneau, Konstantin Kapustin, Randall Leeds, Berker Peksag, Jason Madden, and Brett Randall.
904 gunicorn, 2010. URL <https://gunicorn.org/>. Last accessed 10 July 2022.
- 905 [41] Nginx, Inc. Nginx, 2004. URL <https://nginx.org/en/>. Last accessed 10 July 2022.
- 906 [42] Matthew Bodenhamer. *Learning Relational Concepts with the Spatiotemporal Multidimensional Rela-*
907 *tional Framework*. PhD thesis, University of Oklahoma, 2014.
- 908 [43] Thomas J. Palmer. *Learning Action-State Representation Forests for Implicitly Relational Worlds*. PhD
909 thesis, University of Oklahoma, 2015.
- 910 [44] Peter Pidcoe and Hlapang Kolobe. Self-initiated prone progressive crawler, 01 2015.
- 911 [45] Laura A. Prosser, Julie Skorup, Samuel R. Pierce, Abbas F. Jawad, Andrew H. Fagg, Thubi H. A. Kolobe,
912 and Beth A. Smith. Locomotor Learning in Infants at High Risk for Cerebral Palsy: A Study Protocol.
913 *Frontiers in Pediatrics*, 11, 2023. ISSN 2296-2360. doi: <https://doi.org/10.3389/fped.2023.891633>. URL
914 <https://www.frontiersin.org/articles/10.3389/fped.2023.891633>.
- 915 [46] M Bosanquet, L Copeland, R Ware, and R Boyd. A Systematic Review of Tests to Predict Cerebral
916 Palsy in Young Children. *Developmental Medicine and Child Neurology*, 55(5):418–426, 2013. URL
917 <https://doi.org/10.1111/dmcn.12140>.
- 918 [47] S. K. Campbell. *The Test of Infant Motor Performance. Test User’s Manual Version 3.0 for the TIMP*
919 *Version 5*. Infant Motor Performance Scales, LLC, Chicago, IL, 2012.